


# Roll Adjustable Couch Extension Adaptor



BM02-T manual roll adjustment couch extension adaptor.



BM02-T with SB02 couch extension.



The BM02-T combined with the SB02 couch extension provide both pitch and roll adjustment for head, neck and shoulder patient setup

The BM02-T manual roll adjustment system for Exact style treatment extensions allows for manual adjustment of roll for couch extension accessories. When combined with a CDR SB02 carbon fiber couch extension the combination adds both pitch and roll adjustment to the other axis of a treatment couch.

## BM02-T