




# FREEDOMX™

THE FREEDOMX™ TOTAL BODY PRECISION PATIENT IMMOBILIZATION SYSTEM FOR SRS, SBRT AND MORE.



## A solid foundation is everything.

The FreedomX™ Comprehensive Overlay is the foundation of the FreedomX™ System and has been engineered with you in mind. At 1.7kg/3.8 lbs, this low density overlay weighs less than a bag of sugar but is as strong as steel.


Compatible with all treatment couches using the corresponding indexing bar, the FreedomX™ Comprehensive Overlay is a true precision centric patient positioning system.



Unique 5-point tension lock mechanism ensures masks are securely fastened to the overlay with easy removal.



Slots allow for attachment of FreedomX Intuition Modules (for SRS and Head & Neck) and FreedomX Shoulder Retractors.



Four Point locking system for SBRT components/modules.




Indexing Bars allow indexing to any couch.



Integrated Table Pad ensures patient comfort. Remove and the X feature will accurately index a posterior VacuumLock™.

Complete compatibility with all robotic couch systems including Elekta™ iGuide Tracking System.



<b>FDM-S-STD-X</b>	<b>FreedomX™ System Comprehensive Overlay Package</b> Includes: LB-SBRT-X, SBRT-FILL-X, CFJ02, couch specific carbon fiber index bars (2)	1 pkg.
<b>LB-SBRT-X</b>	FreedomX™ Comprehensive CT/Treatment Overlay, includes CFJ02	1 ea.
<b>SBRT-FILL-X</b>	FreedomX™ Integrated Table Pads	1 ea.
<b>LB-SBRT-WM-X</b>	FreedomX™ Comprehensive Wall Mount	1 ea.
<b>CFJ02</b>	5-Point Tension Lock Immobilization Clamp	1 ea.
<b>RS02</b>	Shoulder Retractors	1 pair
<b>IB06*</b>	Radio-translucent Indexing Bars	2 ea.

*\*Note: contact your sales representative to determine the correct part code for indexing bars that match your couch.*

See [cdrsys.ca](http://cdrsys.ca) for more information on the FreedomX Modules and FreedomX Accessories.