




# FREEDOMX™

THE FREEDOMX™ SBRT MODULE




Thoughtful design and engineering allow you to setup with precision and ease.  
Designed with you in mind!

## Full Body Positioning

As well as Cranial SRS or Head & Neck treatments, the FreedomX™ System can be configured for SBRT type treatment using several additional components that attach and index in conjunction with the FreedomX™ Comprehensive Overlay.

Designed to utilize Vacuum Lock™ cushions, the FreedomX™ System can be manipulated to ensure your patient is extremely comfortable while giving you the highest levels of precision in SBRT patient positioning. The FreedomX™ System is easily configured for Arms Up or Arms Down ensuring your patient will be as comfortable as possible.



The SBRT Module is one of the ten highly-designed FreedomX™ Modules.

**FDM-S-SBRT-X FreedomX™ SBRT Module for Comprehensive Overlay (LB-SBRT-X)** 1pkg.  
Includes: AU-3, KN2-WDG-2, F2-WDG, SBRT-BS-X, SBRT-BL-X, CMP-S-X, CMP-L-X, IB06 type Indexing Bars x2 (specify couch)

**AU-3-SA** FreedomX™ Standalone Arms Up module, indexes directly to most couches 1 ea.  
Contact your sales representative for couch compatibility.

See [cdrsys.ca](http://cdrsys.ca) for more information on the FreedomX Modules and FreedomX Accessories.